

Troubleshooting

Observation	Possible cause	Action	Youtube	Acura Manual
Tips not fitting the nozzle	Inadequat position of tip ejector	Adjust position of tip ejector	https://www.youtube.com/watch?v=ztF42bwEkX0	
	Inadequat tips	use original or other compatible tips		Apply new tip as often as deemed necessary in regard with change of sample nature, droplet accumulation after blow out, or foam formation. When using tip racks, place the instrument on tips, press firmly to ensure a good tightness. This is particularly important for multichannel models. Adequat tips: Qualitix® Pipettenspitze Tips with filter: - Sterile natural ultra-microtip (10mL) - Sterile natural ultra-microtip, extended (20mL) - Sterile natural tip (20mL) - Sterile natural tip (1000mL) - Sterile natural tip (200mL) - Sterile natural tip, extended (200mL) - Sterile natural tip (300mL) - Sterile natural tip (1000 mL) Tips without filter: Vol. (10mL) - Natural ultra-microtip, Natural ultra-microtip, tipfill rack filling system - Empty rack for tipfill system - Natural ultra-microtip, sterile - Natural ultra-microtip, low protein binding - Natural ultra-microtip, low protein binding, sterile Vol. (20mL) - Natural extended microtip - Natural extended microtip, low protein binding - Natural extended microtip, low protein binding, sterile - Natural gel load tip Vol. (10mL) Vol. (200mL) - Natural universal tip - Natural universal tip, tipfill rack filling system - Empty rack for tipfill system - Natural universal tip, sterile - Natural universal tip, low protein binding - Natural universal tip, low protein binding, sterile - Yellow universal tip - Yellow universal tip, sterile - Natural gel load tip - Natural tip Vol. (350mL) Vol. (1000mL) - Natural universal tip - Natural universal tip, tipfill rack filling system - Empty rack for tipfill system - Natural universal tip, sterile - Natural universal tip, low protein binding - Natural universal tip, low protein binding, sterile - Blue universal tip - Blue universal tip, sterile - Natural tip - Blue saver tip - Natural saver tip - Natural macrotip Vol. (2mL) - Natural macrotip Vol. (5mL) - Natural macrotip Vol. (10mL)
Pasteur pipette does not hold or too tight	Pasteur not compatible	Check/change O-rings	https://www.youtube.com/watch?v=0aNLFvVzZpQ&list=PLniRBpWRno86bvIjIoAwHywdUPimX0AqV	Removing volumetric module (lower assembly): - Press ejector button (2), turn part (20) or (24) for 2 mL to the left and pull. - Unscrew barrel (18) - Press plunger button in full then unclip (pull) plunger - After removing volumetric module, take barrel (18), press both clips of ring (14) with a pipette tip or a pointed tool. - Remove ring from barrel, then remove spring (15) to access O-ring/PTFE sleeve assembly. - Change parts if needed. Clean plunger (13) to eliminate any grease residue. Note: To prevent damages when reassembling, first place PTFE sleeve (16) on plunger, then O-ring (17) on sleeve. - Slightly grease plunger, sleeve and O-ring and reassemble volumetric module.
		Grease slightly O-rings		
Pluger sticking or dragging	Dirty plunger	Disassemble and clean plunger	N/A	N/A
Liquid not aspirated	Dry residues block pipette nozzle	Disassemble and clean	N/A	N/A

	Volumetric module not properly assembled	Reassemble according to instructions	N/A	N/A
	Dirty or wet nozzle filter on 835 models	Change nozzle filter	N/A	N/A
Uneven liquid level in tips of multichannel model	Lack of tightness	Check tip fitting	N/A	N/A
		Change for damaged nozzle or tightness parts	N/A	N/A
Poor instrument performance	Lack of tightness	Check for damaged nozzle or tightness parts	N/A	N/A
	Inadequat tips	Check tip compatibility and fitting		Adequat tips: Qualitix® Pipettenspitze Tips with filter: - Sterile natural ultra-microtip (10mL) - Sterile natural ultra-microtip, extended (20mL) - Sterile natural tip (20mL) - Sterile natural tip (1000mL) - Sterile natural tip (200mL) - Sterile natural tip, extended (200mL) - Sterile natural tip (300mL) - Sterile natural tip (1000 mL) Tips without filter: Vol. (10mL) - Natural ultra-microtip, Natural ultra-microtip, tipfill rack filling system - Empty rack for tipfill system - Natural ultra-microtip, sterile - Natural ultra-microtip, low protein binding - Natural ultra-microtip, low protein binding, sterile Vol. (20mL) - Natural extended microtip - Natural extended microtip, low protein binding - Natural extended microtip, low protein binding, sterile - Natural gel load tip Vol. (10mL) Vol. (200mL) - Natural universal tip - Natural universal tip, tipfill rack filling system - Empty rack for tipfill system - Natural universal tip, sterile - Natural universal tip, low protein binding - Natural universal tip, low protein binding, sterile - Yellow universal tip - Yellow universal tip, sterile - Natural gel load tip - Natural tip Vol. (350mL) Vol. (1000mL) - Natural universal tip - Natural universal tip, tipfill rack filling system - Empty rack for tipfill system - Natural universal tip, sterile - Natural universal tip, low protein binding - Natural universal tip, low protein binding, sterile - Blue universal tip - Blue universal tip, sterile - Natural tip - Blue saver tip - Natural saver tip - Natural macrotip Vol. (2mL) - Natural macrotip Vol. (5mL) - Natural macrotip Vol. (10mL)
	Dirty or wet nozzle filter on 835 models	Change nozzle filter	N/A	N/A
	Instrument out of calibration	Perform new calibration	https://www.youtube.com/watch?v=096aGez7Cxc	The swift-set calibration system* with integrated key and locking mechanism makes tedious procedures obsolete. If performance results are no longer within recommended values, for instance after QC check or replacing parts, or if changes occur in physical parameters (liquid density, temperature), recalibration is easily performed using an analytical balance. Ten measurements should be first performed on small volume after checking seal tightness and perfect working condition of the instrument. OR Choice of Service Levels (https://socorex.com/en/laboratory/service-calibration/programme-iso-8655-17025#service-programme-1)

	Instrument used with viscous or volatile liquids or at particular temperature	Perform new calibration with liquid or at temperature used		The swift-set calibration system* with integrated key and locking mechanism makes tedious procedures obsolete. If performance results are no longer within recommended values, for instance after QC check or replacing parts, or if changes occur in physical parameters (liquid density, temperature), recalibration is easily performed using an analytical balance. Ten measurements should be first performed on small volume after checking seal tightness and perfect working condition of the instrument. OR Choice of Service Levels (https://socorex.com/en/laboratory/service-calibration/programme-iso-8655-17025#service-programme-1)
--	---	--	--	--